


Zementrestriktor Cement Restrictor

• Universal-Femur PE-Restriktor / Universal-Femur PE-Restriktor



Zementrestriktor Cement Restrictor



**Zementrestriktor
Cement Restrictor**


Art.-Nr./Art.-No. HT 3500-25

**Applikator für Zementrestriktor
Applicator for Cement Restrictor**

Art.-Nr./Art.-No. IH 0003-00

Der Universal-Femur-Zementrestriktor aus Polyethylen verhindert das Eindringen des Knochenzementes in den distalen Markraum.

The Universal-Femur-Cement Restrictor made out of polyethylene prevents a penetration of bone cement into the distal medullary cavity.



▶ aap Implantate AG
Lorenzweg 5 • 12099 Berlin
Germany
Fon +49 30 75019-0
Fax +49 30 75019-111
customer.service@aap.de
www.aap.de